



Molly®
1.0 mL
Preconfigured 2-Step Solution



Volume 
Viscosity 

Device Type	Disposable autoinjector
Primary Container	1 mL long syringe with staked needle
Dosage Type	Fixed Full contents
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle cover
Injection Feedback	Audible click at beginning and end of injection, viewing window

Handling Instructions



Remove cap



Press and hold against injection site to activate



Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications



The SHL Medical logo and Molly are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, and other jurisdictions. Final product specifications are subject to change based on customer requirements. Contact info@shl-medical.com for more information.

www.shl-medical.com