



Molly® 2.25 mL
Preconfigured 2-Step Solution



Volume 

Viscosity 

Device Type

Disposable autoinjector

Primary Container

2.25 mL pre-filled syringe

Dosage Type

Fixed | Full contents

Route of Administration

Subcutaneous

Needle Insertion

Manual

Activation Method

Needle cover

Injection Feedback

Audible click at beginning and end of injection, viewing window

Handling Instructions



Remove cap



Press and hold against injection site to activate



Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual, and tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic cap removal
- ✓ Robust anti-roll design



The SHL Medical logo and Molly are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, and other jurisdictions. Final product specifications are subject to change based on customer requirements. Contact info@shl-medical.com for more information.

www.shl-medical.com